

LEVEL UP TOES LEVELERS

The relationship between pointe shoes and ballet is in the Só Dança DNA. It is with great pleasure that we present you another product to make your life easier in dance.

Toes Levelers protect your toes where they need it most adapting your feet to fit in all Só Dança pointe shoes styles.

Notice that in **Picture 1** not all the toes touch the pointe shoe bottom, but in **Picture 2** they are leveled, as so the support is homogeneous and the body weight is evenly distributed over all the toes, reducing the pressure on the larger toe.

This makes all the difference to the dancer, because the Levelers favor impact cushioning and alignment of the feet when on pointe, helping to improve the health of the feet and helping to prevent callus formation, abrasions and possible deformation.





Figura 1

Figura 2

WHY TO USE THE SÓ DANÇA TOES LEVELERS?

The proper use of the Levelers favors impact cushioning and feet alignment when you are on pointes helping to keep your feet healthy, helping to prevent calluses, abrasions and possible deformation.

SUGGESTIONS ON HOW TO USE



A LEVELER FOR EACH TOE

There are 5 sizes of Levelers, one for each toe. Each Leveler comes with an extender that can be inserted into the Leveler to increase toe length.

Each pack contains 10 levelers and 10 extenders.

HOW TO USE THE SÓ DANÇA LEVELERS?

Note how your toes are and use the Levelers to try to make them more tapered, as in SUGGESTIONS ON HOW TO USE.

To know which toes should receive the Levelers, it is important to do some testing. Begin by placing a Leveler on the smallest toe next to your largest one.

If necessary, insert the extender to extend the length of the leveler. Then put on the shoe and roll up on pointes. If the pressure on the larger toe is still high, place a Leveler on your next smaller toe, and so on. Repeat the tests until you feel comfortable.

FITTING DIFFERENT TOES SIZES

The elastic nature of the Levelers allows it to adapt to all toes sizes. If necessary, use scissors to reduce the leveler and make a single cut.

If your toes are too short and thin, use the leveler indicated for smaller toes.

LEVELERS PRESERVATION

Keep away from perforating materials and after each use baby powder to maintain the elasticity and soft touch.